


Webcam Eye Tracker Instruction

Note: Identifiable information is currently obscured with a gray box for double-blind purposes.


Section A - Before Task

1. Locate the webcam directory on Desktop and find the app.py file in it.



2. Run app.py file. There may be two ways (Terminal / IDE):

- a. Run app.py using Terminal




```
Windows PowerShell
Copyright (C) Microsoft Corporation. All rights reserved.

Install the latest PowerShell for new features and improvements! https://aka.ms/PSWindows

PS C:\Users\...Data/Local/Programs/Python/Python38/python.exe c:\Users\...
/Desktop/webcam/webcam/app.py
Order of stimuli:[6, 7, 11, 10, 17, 15, 20, 24, 22, 3, 21, 19, 0, 13, 4, 5, 23, 2, 18, 14, 1, 8, 12,
9, 16, 100, 25, 36, 34, 26, 29, 30, 32, 31, 27, 28, 33, 35]
cannot find eyetracker
* Serving Flask app 'app'
* Debug mode: on
WARNING: This is a development server. Do not use it in a production deployment. Use a production WSG
I server instead.
* Running on https://127.0.0.1:9999
Press CTRL+C to quit
* Restarting with stat
Order of stimuli:[21, 7, 17, 22, 23, 6, 12, 24, 18, 14, 9, 3, 0, 13, 4, 20, 11, 10, 15, 5, 2, 8, 19,
16, 1, 100, 25, 32, 33, 26, 28, 36, 31, 30, 29, 35, 34, 27]
cannot find eyetracker
* Debugger is active!
* Debugger PIN: 108-334-022
```


b. Run `app.py` using IDE (e.g., VS Code): you may click the start button on the right top



```
File Edit Selection ... Search
app.py x
C:\X... app.py > ...
30 shuffled.insert(12, 13)
31 shuffled.insert(24, 25)
32 shuffled.insert(12, 0)
33 shuffled.insert(25, 100)
34 # ORDER OF STIMULI #
35 # note: want 10 & 11 to stay consecutive (they use the same tree)
36 if shuffled.index(10) < shuffled.index(11):
37     temp = shuffled[shuffled.index(10) + 1]
38     shuffled[shuffled.index(11)] = temp
39     shuffled[shuffled.index(10) + 1] = 11
40 else:
41     temp = shuffled[shuffled.index(11) + 1]
42     shuffled[shuffled.index(10)] = temp
43     shuffled[shuffled.index(11) + 1] = 10
44
45 print(f'Order of stimuli:{shuffled}')
46
47
48 #
49 # setting up eye-tracker and making sure we can get data from it
50 def eye_tracker_setup():
51     found_eyetrackers = tr.find_all_eyetrackers()
```

3. When you see the following information in your Terminal or IDE Terminal, it means the Webcam Eye Tracker is running successfully:

In Terminal:



```
Windows PowerShell
Copyright (C) Microsoft Corporation. All rights reserved.


Install the latest PowerShell for new features and improvements! https://aka.ms/PSWindows

PS C:\Users\... Data\Local\Programs\Python\Python38\python.exe c:/.../
/Desktop/webcam/webcam/app.py
Order of stimuli:[6, 7, 11, 10, 17, 15, 20, 24, 22, 3, 21, 19, 0, 13, 4, 5, 23, 2, 18, 14, 1, 8, 12,
9, 16, 100, 25, 36, 34, 26, 29, 30, 32, 31, 27, 28, 33, 35]
cannot find eyetracker
* Serving Flask app 'app'
* Debug mode: on
WARNING: This is a development server. Do not use it in a production deployment. Use a production WSG
I server instead.
* Running on https://127.0.0.1:9999
Press CTRL+C to quit
* Restarting with stat
Order of stimuli:[21, 7, 17, 22, 23, 6, 12, 24, 18, 14, 9, 3, 0, 13, 4, 20, 11, 10, 15, 5, 2, 8, 19,
16, 1, 100, 25, 32, 33, 26, 28, 36, 31, 30, 29, 35, 34, 27]
cannot find eyetracker
* Debugger is active!
* Debugger PIN: 108-334-022
```


The screenshot shows a VS Code editor with a file named `app.py` open. The code is a Flask application that sets up a web server on `127.0.0.1:9999`. It creates several directories under `C:\Users\...\\webcam\\webcam_participants\\{pid}`. The terminal at the bottom shows the output of running the application, including the message `* Running on https://127.0.0.1:9999` and a list of stimuli. A red box highlights the URL in the terminal, and a red arrow points to it with the text `Press ctrl and click this link`.

5. Make sure to MAXIMISE the browser window before the next step!

6. You may see this page first. Follow the steps below:




7. Next, you will see the welcome page and a small window on the top left of the page in some seconds, and then you are ready to start:




8. Next, you will be directed to the calibration page and do a calibration yourself.


- How to calibrate yourself: please click all the red circles on the page until they turn yellow. When you complete the calibration, you will receive a notification and see your calibration accuracy on the left top of the page. If the accuracy is **below 70%, please recalibrate!**




- The calibration page looks like this:



- c. After you click all of the 9 circles and make them turn **yellow**, you will see your accuracy score on the left top of the page:




- d. If the score exceeds 70%, you can click the Start button on the top of the page to start your tasks.



- e. You will have **a break** during the task, and you will be directed to the calibration page after the break. So please repeat these procedures when you do the calibration yourself, thank you:)

Section B - After Task

9. After you complete all of the 20 tasks, you will see this completion page and you can then close your browser. Don't forget to **stop** the `app.py` file in the terminal using `ctrl+c`



Thank you so much for your participation again. We really appreciate your time! Please let me know if you have any questions!